



Product Name:

Surgical Face Mask with ear loop

Descriptions:

The Surgical Face mask with ear loop are manufactured and packed in cleanroom environment. Manufactured in ISO factories with in-house laboratory.

Idea for use in medical / pharmaceutical industries / food processing, laboratory and so on.

Features:

High filtration 3-layers fabric

Inner & outer fluid resistant non-woven layers, easy breathable

Extra pliable nose piece ensures snugger fit & reduces fog up of glasses.

No fuzzy tickle feeling on nose & face.

High Filtration Capacity Latex free, fiberglass free With Ties to adjust as you feel more comfortable.

Low resistance to respiration Hypoallergenic.

Specifications:

Material: Non-woven,

Nose Bar: 100%recycling metal bar, 8.5cm

Bacteria Filtration efficiency:97%

Skin irritation: almost negative

Particles: less than 100ea

Weight: 3/g

Type: Ear Loop

Size: 175mm*95mm(length x width)

Color: White/Blue

Packaging:

50pcs per Box

40 boxes per carton

2000pcs pe carton